





Gamma Irradiated Tubing



Reduce the risk of microbial contamination

- Save time—many sizes of packaged and pretreated tubing available for quick shipment!
- Tubing is double bagged prior to gamma irradiation
- Irradiation at 32.5 kGy offers optimal microbial control for the tubing's critical fluid pathway
- Gamma-sensitive labeling confirms treatment
- Select from popular biocompatible tubing formulations—contact us for additional options

Tubing	Formulation	Certifications	ID (Size)	Pack size	Catalog number
	Masterflex® platinum-cured silicone	USP Class V EP FDA REACH RoHS ADCF	1/8" (L/S® 16)	25 ft (7.6 m)/pk	MK-96119-16
			3/16" (L/S® 15)		MK-96119-15
			1/4" (L/S® 24)		MK-96119-24
			3/8" (I/P® 73)		MK-96119-73
			1/2" (I/P® 82)		MK-96119-82
	Masterflex® PharMed® BPT	USP Class VI EP FDA NSF REACH RoHS	1/8" (L/S® 16)	25 ft (7.6 m)/pk	MK-96122-16
			3/16" (L/S® 15)		MK-96122-15
			1/4" (L/S® 24)		MK-96122-24
			3/8" (I/P® 73)		MK-96122-73
			1/2" (I/P® 82)		MK-96122-82
	Masterflex® C-Flex® (Formulation 001)	USP Class VI EP REACH RoHS ADCF	1/8" (L/S® 16)	25 ft (7.6 m)/pk	MK-96123-16
			3/16" (L/S® 15)		MK-96123-15
			1/4" (L/S® 24)		MK-96123-24
			3/8" (I/P® 73)		MK-96123-73
			1/2" (I/P® 82)		MK-96123-82
	C-Flex® Clear (Formulation R70-082)	USP Class VI ADCF GMP	1/8"	25 ft (7.6 m)/pk	MK-96124-05
			1/4"		MK-96124-11
			3/8"		MK-96124-16
			1/2"		MK-96124-18